BMC Infectious Diseases



Open Access Technical advance

Evaluation of a novel real-time PCR test based on the ssrA gene for the identification of group B streptococci in vaginal swabs

Martina Wernecke¹, Ciara Mullen¹, Vimla Sharma³, John Morrison³, Thomas Barry^{1,2}, Majella Maher*¹ and Terry Smith¹

Address: 1 Molecular Diagnostics Research Group, National Centre for Biomedical Engineering Science, Galway, Ireland, 2 Microbiology, School of Natural Sciences, National University of Ireland, Galway, Ireland and ³Department of Obstetrics and Gynaecology, University College Hospital Galway, Ireland

Email: Martina Wernecke - martina.wernecke@nuigalway.ie; Ciara Mullen - ciara.mullen@nuigalway.ie; Vimla Sharma - vimla.sharma@nuigalway.ie; John Morrison - john.morrison@nuigalway.ie; Thomas Barry - thomas.barry@nuigalway.ie; Majella Maher* - majella.maher@nuigalway.ie; Terry Smith - terry.smith@nuigalway.ie

* Corresponding author

Published: 4 September 2009

BMC Infectious Diseases 2009, 9:148 doi:10.1186/1471-2334-9-148

This article is available from: http://www.biomedcentral.com/1471-2334/9/148

This is an Open Access article distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/2.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

Received: 11 February 2009 Accepted: 4 September 2009

© 2009 Wernecke et al; licensee BioMed Central Ltd.

Abstract

Background: Despite the implementation of prevention guidelines, early-onset group B streptococci (GBS) disease remains a cause of neonatal morbidity and mortality worldwide. Strategies to identify women who are at risk of transmitting GBS to their infant and the administration of intrapartum antibiotics have greatly reduced the incidence of neonatal GBS disease. However, there is a requirement for a rapid diagnostic test for GBS that can be carried out in a labour ward setting especially for women whose GBS colonisation status is unknown at the time of delivery. We report the design and evaluation of a real-time PCR test (RiboSEQ GBS test) for the identification of GBS in vaginal swabs from pregnant women.

Methods: The qualitative real-time PCR RiboSEQ GBS test was designed based on the bacterial ssrA gene and incorporates a competitive internal standard control. The analytical sensitivity of the test was established using crude lysate extracted from serial dilutions of overnight GBS culture using the IDI Lysis kit. Specificity studies were performed using DNA prepared from a panel of GBS strains, related streptococci and other species found in the genital tract environment. The RiboSEQ GBS test was evaluated on 159 vaginal swabs from pregnant women and compared with the GeneOhm[™] StrepB Assay and culture for the identification of GBS.

Results: The RiboSEQ GBS test is specific and has an analytical sensitivity of I-10 cell equivalents. The RiboSEQ GBS test was 96.4% sensitive and 95.8% specific compared to "gold standard" culture for the identification of GBS in vaginal swabs from pregnant women. In this study, the RiboSEQ GBS test performed slightly better than the commercial BD GeneOhm™ StrepB Assay which gave a sensitivity of 94.6% and a specificity of 89.6% compared to culture.

Conclusion: The RiboSEQ GBS test is a valuable method for the rapid, sensitive and specific detection of GBS in pregnant women. This study also validates the ssrA gene as a suitable and versatile target for nucleic acid-based diagnostic tests for bacterial pathogens.

Background

Group B Streptococcus (GBS) (Streptococcus agalactiae) is one of the leading causes of neonatal morbidity and mortality in the developed world. Early-onset GBS disease occurs within the first week of life and is associated with neonatal sepsis, pneumonia and meningitis. The mortality rate averages at 6.5% for early-onset GBS cases and for infected preterm infants it rises to 22.7% [1].

Approximately 10 - 40% of pregnant women carry GBS asymptomatically in their vagina or rectum [2-4]. Transmission to the infant occurs vertically during labour via fetal aspiration of infected amniotic fluid or during passage through the birth canal. Since the implementation of intrapartum antibiotic prophylaxis (IAP) for the prevention of GBS, the incidence of early-onset GBS disease in neonates has decreased significantly although resistance to certain antimicrobials has been reported [5]. Studies in the US show that incidence has declined from 1.5/1000 live births in 1990 [6] to 0.32/1000 live births in 2003 [1]. In the UK, early-onset GBS disease is reported to occur in 0.5/1000 births [7].

Practices recommended by public health authorities for the identification of at-risk women who should receive antibiotic treatment during labour vary internationally. Centres for Disease Control and Prevention (CDC) guidelines for the prevention of perinatal GBS disease recommend a universal, microbiological culture-based prenatal screening strategy at 35 - 37 weeks' gestation combined with IAP for GBS-colonised women [2]. Conversely, the risk factor approach relies on the presence of one or more of the following intrapartum factors as indicators for increased risk of neonatal GBS infection: previous infant with early-onset GBS disease, premature labour, prolonged rupture of membranes and fever [7]. A comparison of the two strategies showed that the culture-based approach was almost 50% more effective in preventing early-onset GBS disease [8]. The UK-based Royal College of Obstetricians and Gynaecologists (RCOG) opposes the practice of universal screening and widespread use of intrapartum antibiotics due to a lack of clear evidence of the effectiveness of such practices in controlling the incidence of neonatal sepsis [7]. A recent UK study which looked at the cost-effectiveness of prenatal screening and treatment strategies for the prevention of GBS and other bacterial infections in early infancy proposed that treatment of all pre-term and high-risk term groups would be beneficial [9].

Perinatal GBS disease prevention practices have been successful in lowering neonatal GBS disease by 50 - 80% [10,11]. Nevertheless, early-onset GBS disease cases continue to occur, leading to acute clinical complications especially for preterm infants [12]. There is a requirement

for a rapid diagnostic test for GBS that can be performed in a labour ward setting to ascertain the GBS colonisation status of women in labour, those in preterm labour or women who have not had prenatal care. For these women, culture-screening is not useful because of the time it takes to obtain test results. Several reports demonstrate that real-time PCR is a rapid, more sensitive method than standard culture for determining the intrapartum GBS colonisation status [13-17]. As a result, intrapartum antibiotic prophylaxis can be administered more effectively, thereby reducing the transmission rates of GBS to infants and consequently lowering infant morbidity and mortality rates.

The objective of this study was to develop a rapid nucleic acid diagnostic test for GBS that can be performed in a wide variety of labour and delivery settings. The real-time PCR test developed in this study targets the bacterial *ssrA* gene which was previously shown to be a versatile diagnostic target for important food and clinical pathogens [18,19]. The performance of the *RiboSEQ* GBS test was compared to the commercial FDA-cleared BD GeneOhm™ StrepB Assay (Becton, Dickinson and Company, New Jersey, USA). The results from these tests were correlated with the results of microbiological analyses of the samples.

Methods

Collection of specimens

Vaginal swabs from pregnant women (n=39) were sourced from the Dept. of Obstetrics and Gynaecology, University College Hospital Galway (UCHG), Ireland. Ethics consent was obtained for the study from the Research Ethics Committee at UCHG for vaginal swabs only. Duplicate vaginal swabs were collected into Amies transport medium (Sarstedt, Nümbrecht, Germany), transported to the laboratory at ambient temperature within 3 to 4 hours of collection and processed immediately.

Vaginal swab specimens from pregnant women (n = 120) were purchased from The New England Life Science Group (NELSG) (Los Osos, CA, USA), a clinical services organization. These specimens were remnant swabs screened for GBS colonisation by USA hospital laboratories as part of routine prenatal care. GBS was identified in these swabs at source by genital screen cultures or by selective GBS culture. The remnant specimens were frozen within 1 - 3 days of sampling and shipped on dry ice to our laboratory. A proportion of 90% GBS-positive specimens (n = 107) was requested from NELSG. Those undertaking the performance evaluation of the nucleic acid tests were blinded to the results of the US microbiological analysis while the study was in progress.

Culture

One of the duplicate antenatal specimens collected from UCHG was used for microbiological evaluation. The swab was inoculated into 7 ml LIM broth (Todd-Hewitt broth with 15 μg/ml nalidixic acid and 10 μg/ml colistin) (LIP Diagnostic Services, Galway, Ireland), incubated overnight and subcultured onto TSA + 5% sheep blood agar (LIP) for 18 - 24 h at 37 °C according to CDC recommendations for GBS culture processing. Presumptive GBS colonies showing β-haemolysis were confirmed using the catalase test and antigen detection (Streptococcal Grouping kit, Oxoid, Cambridge, UK). Plates showing no colony growth were reincubated for a further 24 hours and reinspected. Remnant specimens collected from US sites were screened at source for GBS using the CDC-recommended method described above or by genital culture including LIM broth culture.

Sample preparation

The commercially available crude lysis kit BD GeneOhm™ Lysis Kit (BD, NJ, USA) was used for sample preparation. The second swab of duplicate UCHG specimens and remnant swabs were vortexed for 2 min in 1 ml of sample buffer. Four hundred microlitres of the suspension was transferred into the lysis tube and lysed by mechanical disruption with silica beads according to manufacturer's instructions. Lysates were stored at -20°C until required.

Construction of IAC

An internal amplification control (IAC) was constructed using the composite primer approach described by Hoorfar et al. [20,21] and O'Grady et al. [19]. An internal amplification control consisting of heterologous DNA cloned into a plasmid vector was included in the *RiboSEQ* GBS test to identify false negative test results caused by PCR inhibition. The IAC is co-amplified with the GBS target and detected by an IAC-specific hybridization probe (Table 1) included in the real-time PCR reaction.

RiboSEQ GBS real-time PCR test

The ssrA genes of 10 GBS strains which represent 7 serotypes including the 5 most commonly occurring serotypes (Ia, Ib, II, III, V) were sequenced (Sequiserve, Vaterstetten, Germany) and aligned with other related streptococci (Table 2) ssrA sequences generated in this study or available on the tmRNA website [22]. From these alignments, oligonucleotide primers gbsU3F and gbsU4R were designed to amplify a 293 bp PCR product from the GBS ssrA gene. A fluorescently labelled hybridization probe pair (f1GBS-Flu and f2GBS-LC640) was designed for the detection of GBS. BLAST (Basic Local Alignment Search Tool) analysis of the hybridization probe sequences was performed to confirm in silico specificity of the probes for the detection of GBS.

Real-time PCR was performed in a 20 μ l reaction volume on the LightCycler® instrument using the "LightCycler® FastStart DNA Master HybProbe" kit (Roche Diagnostics, Mannheim, Germany). Each reaction contained reagents to final concentrations of: 5 mM MgCl₂, 0.5 μ M each primer (Table 1), 0.2 μ M each hybridization probe (Table 1), 0.5 U uracil-DNA glycosylase (Roche). Template was added in 2 μ l volumes and IAC was added as 100 recombinant plasmid copies per reaction. Thermal cycling parameters consisted of 95 °C denaturation for 10 min followed by 50 amplification cycles of 95 °C for 10 s, 50 °C for 15 s and 72 °C for 10 s. Melting profiles were run between 40 °C and 80 °C at a transition rate of 0.1 °C/s.

Analytical sensitivity of the RiboSEQ GBS real-time PCR test

The limit of detection (LOD) was established using crude lysate extracted from serial dilutions of overnight GBS culture (BCCM 15081) using the IDI Lysis kit (GeneOhm Sciences, Canada). Colony forming units (cfu) per dilution were established by triplicate plate counts and real-time PCR reactions included GBS crude lysate template containing between 10⁵ and 10⁻¹ cell equivalents.

Table I: Oligonucleotide primers and hybridization probes used in this study

Name	Function	Sequence 5' - 3'		
gbsU3F	Forward primer	GACAGGCATTATGAGGTA		
gbsU4R	Reverse primer	GCTAATATTTGTCTACAAC		
CompF	Forward composite primer for IAC generation	GACAGGCATTATGAGGTAATACCCAACTTGGAATG		
CompR	Reverse composite primer for IAC generation	GCTAATATATTTGTCTACAACTCTTCACCAGAATAAAATTG		
tmUF	Universal bacterial ssrA forward primer	GGGG(A/C)(C/T)TACGG(A/T)TTCGAC		
tmUR	Universal bacterial ssrA reverse primer	GGGA(A/G)TCGAACC(A/G)(C/G)GTCC		
f I GBS-Flu	GBS-specific hybridization probe	TTGCGTTTTGCTAGAAGGTCTTA - Flu		
f2GBS-LC640	GBS-specific hybridization probe	LC640- TATCAGCAAACTACGTTTGGCT - Ph		
ALS1-Flu	IAC-specific hybridization probe	TGAATGTATCCCCTGGA - Flu		
ALSI-LC705	IAC-specific hybridization probe	LC705 - TGGCACTGGTACCATCTAA - Ph		

Table 2: Bacterial species and strains used to test the specificity of the RiboSEQ real-time PCR test

GBS	Serotype	Collection ref no	Stool Panel	Collection ref no
Streptococcus agalactiae	la	BCCM 15081	Escherichia coli	DSM 30038
Streptococcus agalactiae	lb	BCCM 15082	Morganella morganii	DSM 30164
Streptococcus agalactiae	lc	BCCM 15083	Proteus mirabilis	DSM 4479
Streptococcus agalactiae	II	BCCM 15084	Klebsiella pneumoniae	Clinical isolate
Streptococcus agalactiae	III	BCCM 15085	Escherichia hermannii	DSM 4560
Streptococcus agalactiae	IV	BCCM 15086	Escherichia vulneris	DSM 4564
Streptococcus agalactiae	٧	BCCM 15087	Aeromonas hydrophila	DSM 30015
Streptococcus agalactiae	lb	BCCM 15090	Citrobacter freundii	DSM 30039
Streptococcus agalactiae	III	BCCM 15094	Enterobacter cloacae	ATCC 13047
Streptococcus agalactiae	III	BCCM 15095		
,			Respiratory Panel	
Strep spp			Staphylococcus aureus	DSM 11965
Streptococcus dysgalactiae		DSM 6176	Staphylococcus epidermidis	ATCC 27626
Streptococcus oralis		DSM 20627	Pseudomonas aeruginosa	ATCC 10145
Streptococcus pneumoniae		DSM 11865	Klebsiella oxytoca	DSM 5175
Streptococcus parasanguinis		DSM 6778	•	
Streptococcus intermedius		DSM 20573	Neisseria Panel	
Streptococcus salivarius		DSM 20560	Neisseria meningitidis	DSM 10036
Streptococcus uberis		DSM 20569	Neisseria lactamica	DSM 4691
Streptococcus mitis		DSM 12643	Neisseria meningitidis	ATCC 13077
Enterococcus faecalis		DSM 20371	Neisseria cinerea DSM 4630	
Enterococcus faecium		DSM 20477		
Streptococcus gordonii		DSM 6777	Sepsis Panel	
Streptococcus mutans		DSM 20523	Streptococcus sanguis	DSM 20567
Streptococcus pyogenes		DSM 20565	Streptococcus porcinus	DSM 20725
Streptococcus anginosus		DSM 20563	Serratia marcescens	DSM 1608
			Enterobacter aerogenes	DSM 30053
Genital Panel			Staphylococcus haemolyticus	DSM 20263
Acinetobacter Iwoffii		DSM 2403	Streptococcus bovis	DSM 20480
Lactobacillus acidophilus		DSM 20079	•	
Lactobacillus fermentum		DSM 20055		
Mobiluncus curtisii		DSM 2711		
Mobiluncus mulieris		DSM 2710		

BD GeneOhm™ StrepB Assay

The performance of the RiboSEQ GBS real-time PCR test for the detection of GBS in clinical samples was benchmarked against the performance of the commercially available FDA-approved BD GeneOhm™ StrepB Assay which targets the CAMP factor B (cfb) gene. The assay was performed on the SmartCycler® instrument (Cepheid, California, USA) according to manufacturer's instructions. Briefly, 25 µl of diluent was added to each Master Mix reaction tube. Positive and negative control tubes were included in each run. Sample crude lysates were added to the Master Mix tubes in 1.5 µl volumes, the tubes were briefly centrifuged and placed in the SmartCycler® instrument. Assay runs were performed using the BD Gene-Ohm™ StrepB Assay-specific software. The GeneOhm™ StrepB Assay takes approximately 60 min to run including sample preparation.

Results

Analytical sensitivity and specificity of the RiboSEQ GBS real-time PCR test

The limit of detection of the *RiboSEQ* GBS test was determined by three independent experimental assessments and consistently detected in the range of 1 - 10 cell equivalents per reaction. The inclusion of an internal amplification control at a concentration of 100 plasmid copies per test did not lower the LOD (data not shown). At high GBS template concentrations ($\geq 10^5$ cell equivalents/reaction) the amplification of internal control target was occasionally inhibited due to competition between GBS target DNA and IAC in the PCR reaction.

Specificity studies were performed using DNA prepared from a panel of 10 GBS strains representing 7 serotypes including 5 common serotypes and 42 related strepto-cocci and other species found in the genital tract environment (Table 2). All 10 GBS strains were detected. The test was negative for all non-GBS species listed in Table 2. The DNA extracted from non-GBS species was confirmed to be

amenable to real-time PCR with universal primers tmUF and tmUR (Table 1) which were designed to amplify the bacterial *ssrA* gene.

Evaluation of the RiboSEQ GBS test in clinical samples

Performance of the *RiboSEQ* GBS test was compared to the results of microbiological culture methods for the isolation of GBS. Test results were also compared to those obtained using the commercial BD GeneOhm $^{\text{TM}}$ StrepB Assay. Microbiological culture results were used as the gold standard.

A total of 159 specimens (39 from UCHG, 120 from NELSG) were tested by microbiological culture, *RiboSEQ* GBS real-time PCR and BD GeneOhm™ StrepB Assay. Microbiological culture identified 111 samples as GBS-positive and 48 samples as GBS-negative.

The prevalence of GBS colonisation among the specimens collected from UCHG was 10.3% (4/39) using the GBS selective culture method. Calculation of prevalence of GBS colonisation in the NELSG samples is not applicable since this population of specimens was selected to be predominantly GBS-positive and not representative of the true prevalence among pregnant women.

In 146 samples, the results of all three tests showed 100% correlation, where 104 samples were identified as GBS-positive and 42 samples were identified as GBS-negative. In 95% of samples the two nucleic acid diagnostic test (NAD) results agreed (151/159).

Discrepant results were obtained in a total of 13 samples after retesting the crude lysate. In 5 samples, NAD test results agreed but conflicted with the microbiology results. Microbiological culture identified 2 samples as GBS-negative which the NAD tests identified as GBS-positive. Three samples were determined as GBS-positive by culture and tested negative in the NAD tests.

In one sample, the *RiboSEQ* GBS test was negative for a sample that was GBS-positive by both culture and BD GeneOhm[™] StrepB Assay. PCR inhibition was not apparent since the IAC gave a positive signal.

For 6 samples, the BD GeneOhm[™] StrepB Assay results disagreed with concurring culture and GBS real-time PCR results (3 false positives, 3 false negatives).

In another sample, both the IAC and the *RiboSEQ* GBS test results were negative, indicating PCR inhibition. This specimen was identified as GBS-negative by both culture and BD GeneOhmTM StrepB Assay with no indication of PCR inhibition. However, the discrepancy was resolved for this sample by increasing the IAC to 200 plasmid copies in the *RiboSEQ* GBS real-time PCR reaction which yielded a positive result for the IAC for this sample.

Table 3 shows the sensitivity, specificity, predictive values and likelihood ratios for the RiboSEQ GBS test and the BD GeneOhmTM StrepB Assay using microbiological culture for the identification of GBS as the gold standard. Confidence intervals (CI) are stated at the 95% level.

Table 3: Sensitivity, specificity, predictive values (PPV, NPV) and likelihood ratios +LR, -LR) for the RiboSEQ GBS test and the BD GeneOhm™ StrepB Assay

		RiboSEQ GBS test		BD GeneOhm [™] StrepB Assay	
		Positive	Negative	Positive	Negative
Culturea	Positive	107	4	105	6
	Negative	2	46	5	43
	Sensitivity	96.4%		94.6%	
	95% CI	(90.5 - 98.8%)		(88.1 - 97.8%)	
	Specificity	95.8%		89.6%	
	95% CI	(84.6 - 99.3%)		(76.6 - 96.1%)	
	PPV ^b	98.2%		95.5%	
	95% CI	(92.9 - 99.7%)		(89.2 - 98.3%)	
	NPVb	92%		87.8%	
	95% CI	(79.9 - 97.4%)		(74.5 - 94.9%)	
	+LR ^b	23.1		9.1	
	95% CI	(5.9 - 89.9)		(3.9 - 20.8)	
	-LR ^b	0.04		0.06	
	95% CI	(0.01 - 0.1)		(0.03 - 0.13)	

Calculated using programme available at http://faculty.vassar.edu/lowry/clin1.html

^a using culture as the gold standard (n = 159)

^b sample population comprised of random (n = 39) and non-random (n = 120) specimens

Discussion

In a systematic review of the accuracy and rapidity of various intrapartum GBS colonisation tests, real-time PCR was identified as one of two feasible methods for maternal intrapartum screening [23]. Several studies have demonstrated the speed, sensitivity and accuracy of real-time PCR tests for the detection of GBS in pregnant women [14,24,25]. Recently, a real-time PCR assay for the detection of ST-17 serotype III, a virulent strain of GBS associated with neonatal invasive infections was developed [26]. Clinical studies have shown an increase in sensitivity for real-time PCR tests in comparison to the gold standard culture method [13,27]. Higher detection rates of realtime PCR tests may be explained by the inability of culture to detect low numbers of organisms, the presence of antagonistic organisms, or possibly the detection of nonviable cells [10].

A relatively low carriage rate obtained by culture for GBS in this study was similar to the carriage rate of 8% reported in a recent study where vaginal swabs from pregnant women were screened for GBS [28]. Ideally, both vaginal and rectal swabs are required to ensure the highest possible GBS detection rate [29]. However, for this study rectal samples were not available since the specimens were collected during routine prenatal care for the purpose of general microbiological screening. Additionally, only colonies showing β -haemolysis were tested with the catalase and latex agglutination tests. The CDC method also recommends testing characteristic colonies without β -haemolysis.

In this study we developed a qualitative real-time PCR test for the detection of GBS that is rapid (75 min including sample preparation) and capable of detecting 1 - 10 genome copies of GBS. The sensitivity and specificity achieved with the test in comparison to culture were 96.4% (CI 95% 90.5 - 98.8) and 95.8% (CI 95% 84.6 -

99.3), respectively. These values compare well with previously published sensitivities and specificities of real-time PCR assays for GBS [14,15,29,24,27] (Table 4). In our hands, the *RiboSEQ* GBS real-time PCR test performed slightly better than the commercial FDA-approved BD GeneOhm™ StrepB Assay (sensitivity 94.6% (CI 95% 88.1 - 97.8), specificity 89.6% (CI 95% 76.6 - 96.1)). Two previous studies evaluated the clinical performance of the BD GeneOhm™ StrepB Assay and reported a sensitivity of 94% and a specificity of 95.9% for direct detection [14] and a sensitivity of 92.5% and a specificity of 92.5% after 4 h selective enrichment [16].

Test results from the real-time PCR tests correlated in 95% of samples. In five samples the culture results disagreed with concurring NAD test results. Two of these 5 samples were identified as culture-negative but GBS-positive by both real-time PCR tests. Increased sensitivities for realtime PCR diagnostic tests over the standard culture method have been reported in other studies [8,13,14] and can be a result of detection of non-viable cells, low bacterial burden [13,16] or the presence of antagonistic microorganisms which inhibit growth in culture [30,31]. In three cases, the PCR tests were negative for culture-positive samples. However, this discrepancy may be due to misidentification of GBS during culture screening. Two of these specimens were screened by genital culture and GBS growth may have been misidentified. The third sample was a UCHG duplicate swab in which the DNA may have been degraded. This sample had been stored at 4°C for 5 days before processing for real-time PCR.

Due to the dynamic status of vaginal GBS colonisation, screening intrapartum is the most accurate method of predicting the GBS colonisation status [2,32]. Screening for GBS at the time of delivery allows antibiotics to be used more effectively thereby rationalising their use and minimising associated risks. The value of a rapid and accurate

Table 4: Comparison of sensitivities and specificities of a selection of diagnostic real-time PCR assays for GBS (using culture for GBS as gold standard)

Reference	Assay name	Target gene	Sensitivity % (95% CI)	Specificity % (95% CI)
This study	RiboSEQ GBS	ssrA	96.4 (90.5 - 98.8)	95.8 (84.6 - 99.3)
Bergeron et al (2000)	Not specified	cfb	97 (82.5 - 99.8)	100 (94.2 - 100)
Bergseng et al (2007)	Not specified	sip	97 (90 - 99)	99 (97 - 100)
Davies et al (2004)	IDI-Strep B	cfb	94 (90.1 - 97.8)	95.9 (94.3 - 97.4)
Gavino and Wang (2007)	GeneXpert GBS	cfb	95.8 (76.9 - 99.8)	64.5 (45.4 - 80.2)
Uhl et al (2005)	LightCycler Strep B	ptsl	100 (88.8 - 100)	96.9 (92.2 - 99.1)

test for GBS colonisation is especially apparent in cases of premature delivery of <35 weeks gestation when screening results may not be available and babies are at greatest risk of developing early-onset GBS disease.

In spite of the advantages of sensitivity, specificity and fast turnaround time, the complexity of real-time PCR technology poses limitations to its widespread application in clinical laboratory testing. Generally, it requires trained personnel and dedicated laboratory areas. In addition, specimens are often batched for analysis prolonging turnaround time. However, these limitations have been reduced through the development of fully integrated and automated *in vitro* diagnostic platforms which are moderately complex in their operation and can be performed in an on-demand schedule [33]. These platforms provide the potential for diagnostic tests such as the *RiboSEQ* GBS test to be employed in hospital laboratories and near-patient settings where a rapid and reliable diagnosis is most valuable.

Conclusion

In summary, a rapid and sensitive qualitative test for the detection of GBS in vaginal swabs by real-time PCR was developed. The test targets the bacterial *ssrA* gene which was previously shown to be a suitable and versatile diagnostic target for important food and clinical pathogens. The method has potential to be employed as a screening and/or diagnostic test for use in clinical laboratories where an assessment of the patient's GBS colonisation status at the time of delivery is required without delay.

Competing interests

Some of the authors are inventors of a platform technology based on the *ssrA* gene (*RiboSEQ*) for bacterial identification that is currently going through PCT patent filing.

Authors' contributions

MW: Performed the experimental work described in the manuscript and contributed intellectually to the design of the real-time PCR RiboSEQ GBS test. Lead author in preparing the manuscript. CM: Performed the experimental work described in the manuscript and contributed to manuscript preparation. VS: Contributed to clinical study design, collection of specimens for the study and analysis and interpretation of results. JM: Contributed to clinical study design, collection of specimens for the study and critical analysis and interpretation of results. TB: Member of the project management team responsible for the study concept and design, critical review of the experimental work and review of the manuscript. MM: Member of the project management team responsible for the study concept and design, critical review of the experimental work and preparation and critical review of the manuscript. TS: Member of the project management team responsible for the study concept and design, critical review of the experimental work and preparation and critical review of the manuscript.

Acknowledgements

This work was supported by the Enterprise Ireland Commercialisation Fund Technology Development Programme (Grant no. CFTD 2005/112).

References

- Centers for Disease Control and Prevention (CDC): Diminishing racial disparities in early-onset neonatal group B streptococcal disease--United States, 2000-2003. Morb Mortal Wkly Rep 2004. 53:502-505.
- Centres for Disease Control and Prevention (CDC): Prevention of perinatal group B streptococcal disease. Morb Mortal Wkly Rep 2002, 51:1-6.
- Meyn LA, Moore DM, Hillier SL, Krohn MA: Association of sexual activity with colonisation and vaginal acquisition of group B Streptococcus in nonpregnant women. Am J Epidemiol 2002, 155:949-957.
- Elbaradie SMY, Mahmoud M, Farid M: Maternal and neonatal screening for group B Streptococci by scpB gene based PCR: A preliminary study. Ind I Med Microbiol 2009. 27:17-21.
- A preliminary study. Ind J Med Microbiol 2009, 27:17-21.

 5. Castor ML, Whitney GC, Como-Sabetti K, Facklam RR, Ferrieri P, Bartkus JM, Juni BA, Cieslak PR, Farley MM, Dumas NB, Schrag SJ, Lynfield R: Antibiotic resistance pattern in invasive group B streptococcal isolates. Infect Dis Obstet Gynecol 2008.
- Zangwill KM, Schuchat A, Wenger JD: Group B streptococcal disease in the United States, 1990: report from a multistate active surveillance system. Morb Mortal Wkly Rep 1992, 41(SS-6):25-32.
- Royal College of Obstetricians and Gynaecologists (RCOG): Prevention of Early-onset Neonatal Group B Streptococcal Disease. In Green-top Guideline Volume 36. London: RCOG; 2003.
- Schrag SJ, Zell ER, Lynfield R, Roome A, Arnold KE, Craig AS, Harrison LH, Reingold A, Stefonek K, Smith G, Gamble M, Schuchat A:
 Active Bacterial Core Surveillance Team. A population-based comparison of strategies to prevent early-onset group B streptococcal disease in neonates. N Engl J Med 2002, 347:233-239.
- Colbourn T, Asseburg C, Bojke L, Philips Z, Claxton K, Ades AE, Gilbert RE: Prenatal screening and treatment strategies to prevent group B streptococcal and other bacterial infections in early infancy; cost effectiveness and expected value of information analysis. Health Technol Assess 2007, 11:1-226iii.
- Daley AJ, Garland SM: Prevention of neonatal Group B streptococcal disease: Progress, challenges and dilemmas. J Paediatr Child Health 2004, 40:664-668.
- Schrag SJ, Zywicki S, Farley MM, Reingold AL, Harrison LH, Lefkowitz LB, Hadler JL, Danila R, Cieslak PR, Schuchat A: Group B streptococcal disease in the era of intrapartum antibiotic prophylaxis. N Engl | Med 2000, 342:15-20.
- Puopolo KM, Madoff LC, Eichenwald EC: Early-Onset Group B Streptococcal Disease in the Era of Maternal Screening. Pediatrics 2005, 115:1240-1246.
- Convert M, Martinetti Lucchini G, Dolina M, Piffaretti JC: Comparison of LightCycler PCR and culture for detection of group B streptococci from vaginal swabs. Clin Microbiol Infect 2005, 11:1022-1026.
- Davies HD, Miller MA, Faro S, Gregson D, Kehl SC, Jordan JA: Multicenter study of a rapid molecular-based assay for the diagnosis of group B Streptococcus colonization in pregnant women. Clin Infect Dis 2004, 39:1129-1135.
- Gavino M, Wang E: A comparison of a new rapid real-time polymerase chain reaction system to traditional culture in determining group B streptococcus colonisation. Am J Obstet Gynecol 2007, 197:e1-388e4.
- Goodrich JS, Miller MB: Comparison of culture and 2 real-time polymerase chain reaction assays to detect group B Streptococcus during antepartum screening. Diagn Microbiol Infect Dis 2007. 59:17-22.
- Rallu F, Barriga P, Scrivo C, Martel-Laferrière V, Laferrière C: Sensitivities of antigen detection and PCR assays greatly

- increased compared to that of the standard culture method for screening for group B Streptococcus carriage in pregnant women. J Clin Microbiol 2006, 44(3):725-728.
- Glynn B, Lahiff S, Wernecke M, Barry T, Smith TJ, Maher M: Current and emerging molecular diagnostic technologies applicable to bacterial food safety. Int J Dairy Technol 2006, 59:126-139.
 O'Grady J, Sedano-Balbás S, Maher M, Smith T, Barry T: Rapid real-
- O'Grady J, Sedano-Balbás S, Maher M, Smith T, Barry T: Rapid realtime PCR detection of Listeria monocytogenes in enriched food samples based on the ssrA gene, a novel diagnostic target. Food Microbiol 2008, 25(1):75-84.
- Hoorfar J, Ahrens P, Rådström P: Automated 5' nuclease PCR assay for identification of Salmonella enterica. J Clin Microbiol 2000, 38(9):3429-3435.
- Hoorfar J, Malorny B, Abdulmawjood A, Cook N, Wagner M, Fach P: Practical considerations in design of internal amplification controls for diagnostic PCR assays. J Clin Microbiol 2004, 42:1863-1868.
- Williams KP: The tmRNA website. Nucleic Acids Res 2000, 28:168-170.
- 23. Honest H, Sharma S, Khan KS: Rapid tests for Group B Streptococcus colonisation in laboring women: a systematic review. *Pediatrics* 2006, 117:1055-1066.
- Bergeron MG, Ke D, Ménard C, Francois FJ, Gagnon M, Bernier M, Ouellette M, Roy PH, Marcoux S, Fraser WD: Rapid detection of group B streptococci in pregnant women at delivery. N Engl J Med 2000, 343:175-179.
- Ke D, Menard C, Picard FJ, Boissinot M, Ouellette M, Roy PH, Bergeron MG: Development of conventional and real-time PCR assays for the rapid detection of group B streptococci. Clin Chem 2000, 46(3):324-331.
- Lamy MC, Dramsi Š, Billoët A, Réglier-Poupet H, Tazi A, Raymond J, Guérin F, Couvé E, Kunst F, Glaser P, Trieu-Cuot P, Poyart C: Rapid detection of the "highly virulent" group B Streptococcus ST-17 clone. Microbes Infect 2006, 8:1714-1722.
- Uhl JR, Vetter EA, Boldt KL, Johnston BW, Ramin KD, Adams MJ, Ferrieri P, Reischl U, Cockerill FR: Use of the Roche LightCycler Strep B Assay for detection of Group B Streptococcus from vaginal and rectal swabs. | Clin Microbiol 2005, 43:4046-4051.
- Van der Mee-Marquet N, Jouannet C, Domelier AS, Arnault L, Lartique MF, Quentin R: Genetic diversity of Streptococcus agalactiae strains and density of vaginal carriage. J Med Microbiol 2009, 58:169-173.
- Bergseng H, Bevanger L, Rygg M, Bergh K: Real-time PCR targeting the sip gene for detection of group B streptcoccus colonisation in pregnant women at delivery. J Med Microbiol 2007, 56:223-228.
- Dunne WM Jr, Holland-Staley CA: Comparison of NNA agar culture and selective broth culture for detection of group B streptococcal colonisation in women. J Clin Microbiol 1998, 36:2298-2300.
- Ostroff RM, Steaffens JW: Effect of specimen storage, antibiotics, and feminine hygiene products on the detection of group B streptococcus by culture and the STREP B OIA test. Diagn Microbiol Infect Dis 1995, 22:253-259.
- 32. Yancey MK, Schuchat A, Brown LK, Ventura VL, Markenson GR: The accuracy of late antenatal screening cultures in predicting genital group B streptococcal colonisation at delivery. Obstet Gynecol 1996, 88:811-815.
- Boschert S: Fast, point-of-care Group B Strep test gets FDA approval: labor and delivery nurses can perform it. Ob Gyn News 2006, 41:1.

Pre-publication history

The pre-publication history for this paper can be accessed here:

http://www.biomedcentral.com/1471-2334/9/148/prepub

Publish with **Bio Med Central** and every scientist can read your work free of charge

"BioMed Central will be the most significant development for disseminating the results of biomedical research in our lifetime."

Sir Paul Nurse, Cancer Research UK

Your research papers will be:

- available free of charge to the entire biomedical community
- peer reviewed and published immediately upon acceptance
- cited in PubMed and archived on PubMed Central
- ullet yours you keep the copyright

Submit your manuscript here: http://www.biomedcentral.com/info/publishing_adv.asp

