

RESEARCH

Open Access



# Short-term culture for rapid identification by mass spectrometry and automated antimicrobial susceptibility testing from positive bottles

Peng-Peng Tian<sup>1</sup>, Shan-Shan Su<sup>1</sup>, Li-Sha Zhu<sup>1</sup>, Tian Wang<sup>1</sup>, Hui Yang<sup>1</sup>, Meng-Yao Du<sup>1</sup>, Cai-Zhi Ding<sup>3</sup>, Li Wang<sup>1</sup>, Wen Fan<sup>1\*</sup> and Hua-Wei Yi<sup>1,2\*</sup>

## Abstract

**Background** Early and appropriate antibiotic treatment improves the clinical outcome of patients with sepsis. There is an urgent need for rapid identification (ID) and antimicrobial susceptibility testing (AST) of bacteria that cause bloodstream infection (BSI). Rapid ID and AST can be achieved by short-term incubation on solid medium of positive blood cultures using MALDI-TOF mass spectrometry (MS) and the BD M50 system. The purpose of this study is to evaluate the performance of rapid method compared to traditional method.

**Methods** A total of 124 mono-microbial samples were collected. Positive blood culture samples were short-term incubated on blood agar plates and chocolate agar plates for 5 ~ 7 h, and the rapid ID and AST were achieved through Zybio EXS2000 MS and BD M50 System, respectively.

**Results** Compared with the traditional 24 h culture for ID, this rapid method can shorten the cultivation time to 5 ~ 7 h. Accurate organism ID was achieved in 90.6% of Gram-positive bacteria (GP), 98.5% of Gram-negative bacteria (GN), and 100% of fungi. The AST resulted in the 98.5% essential agreement (EA) and 97.1% category agreements (CA) in NMIC-413, 99.4% EA and 98.9% CA in PMIC-92, 100% both EA and CA in SMIC-2. Besides, this method can be used for 67.2% (264/393) of culture bottles during routine work. The mean turn-around time (TAT) for obtaining final results by conventional method is approximately  $72.6 \pm 10.5$  h, which is nearly 24 h longer than the rapid method.

**Conclusions** The newly described method is expected to provide faster and reliable ID and AST results, making it an important tool for rapid management of blood cultures (BCs). In addition, this rapid method can be used to process most positive blood cultures, enabling patients to receive rapid and effective treatment.

**Keywords** Bloodstream infection, Short-term culture, Turn-around time, Rapid identification, Antimicrobial susceptibility testing

\*Correspondence:

Wen Fan

931412421@qq.com

Hua-Wei Yi

yihuawei0902@126.com

<sup>1</sup>Laboratory Department, The First Affiliated Hospital of Yangtze University, Jing Zhou, Hubei, China

<sup>2</sup>Hubei Provincial Clinical Research Center for Individualized Diagnosis and Treatment of Cancer, The First People's Hospital of Jingzhou, Jingzhou, Hubei, China

<sup>3</sup>Laboratory Department, The People's Hospital of Songzi, Jingzhou, Hubei, China



© The Author(s) 2024. **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit <http://creativecommons.org/licenses/by/4.0/>. The Creative Commons Public Domain Dedication waiver (<http://creativecommons.org/publicdomain/zero/1.0/>) applies to the data made available in this article, unless otherwise stated in a credit line to the data.

## Introduction

Bloodstream infection (BSI) is one of the leading causes of death worldwide [1, 2]. BSI has caused a significant global healthcare burden, with an estimated mortality rate of 12~20% in 2017 [3]. It is well established that the survival of BSI patients depends on the rapid administration of effective antimicrobial therapy [4, 5]. The survival rate decreases by approximately 7.6% when antimicrobial administration is delayed for 1 h [6]. Rapid identification (ID) and antimicrobial susceptibility testing (AST) of the BSI pathogens are essential for timely selection of appropriate antimicrobial therapy, which may result in a better outcome for patients [7, 8]. Blood culture (BC) is the gold standard method for the diagnosis of BSI, which includes processes such as sample collection, incubation, ID, and AST.

The conventional method requires positive BCs to be cultured on solid medium for 18~24 h or more before ID and AST. Matrix-assisted laser desorption/ionization time-of-flight mass spectrometry (MALDI-TOF MS) has proven to be a rapid, accurate, and cost-effective technology in the routine identification of microorganisms [9, 10]. The mass spectrometer can provide excellent identification in positive blood cultures and has also been applied in drug resistance, such as direct detection of *Klebsiella pneumoniae* carbapenemase (KPC) based on MALDI-TOF MS [11, 12]. However, these methods require additional labor-intensive, expensive steps in ID, and cannot effectively shorten the time of AST.

In order to reduce the turn-around time (TAT) of current methods, Jihye Ha et al. established a rapid method by short-term incubation of positive blood culture samples on solid culture medium for 6 h followed by ID and AST [13]. Here, this study aims to evaluate the performance of the rapid method. After a 5~7 h incubation of positive blood culture samples, MALDI-TOF MS and automated devices were used to detect ID and AST, respectively. In addition, we hope to integrate this method into our routine laboratory to shorten the incubation time for ID and AST, thereby optimizing the standard procedures in microbiology laboratory.

## Methods

### Blood culture samples

Blood culture (BC) bottles (Bactec plus/F; Becton Dickinson, Franklin Lakes, NJ, USA) were collected from patients with suspected BSI between August to September of 2023 at the First Affiliated Hospital of Yangtze University, Jing Zhou, China. The BC bottles were incubated in the Bactec system (Becton Dickinson) until a positive result was obtained or for a maximum of 5 days. In this study, only positive BC samples with Gram staining as distinctly single microorganisms were selected. A total of 124 positive BCs were analyzed using the conventional

laboratory diagnostic method and the novel rapid method. If the positive BCs were processed before 10 am, colonies incubated on solid culture medium for about 5~7 h could undergo rapid ID and AST around 3:30 pm in the afternoon.

### Conventional method of ID and AST

After the BD Bactec system displayed a positive signal, Gram staining was performed, and then sub-cultured on blood agar plates (BAP) and chocolate agar plates (CAP) containing vancomycin. These plates were grown in an incubator (Thermo Fisher Scientific, USA) with 5% CO<sub>2</sub> at 37°C. After overnight incubation, the colonies grown on the BAP were used for ID using MALDI-TOF MS (Zybio EXS2000, China). Bacterial colonies were transferred to the MS target plates using a wooden toothpick. The target plates were overlaid with 1 µL of 70% formic acid. Once the formic acid solution dried, 1 µL of matrix solution was added for subsequent MALDI-TOF MS ID. The calibration and validation of MALDI-TOF MS were carried out once a week with a bacterial test standard according to the manufacturer's instructions. In brief, a score ≥ 2.0 was interpreted as reliable ID to the species level, a score of 1.7-2.0 was interpreted as reliable ID to the genus level, and a score < 1.7 was interpreted as no reliable ID.

A standardized inoculum (McFarland standard of 0.5) was then prepared from single colonies grown on the BAP, and the appropriate BD Phoenix™ M50 AST panels were chosen according to the ID results provided by Zybio EXS2000 MS. AST plates NMIC-413, PMIC-92, and SMIC/ID-2 were used for Gram negative (GN) bacteria, *Staphylococcus/Enterococcus* spp, and *Streptococci* spp, respectively. AST results were obtained after 18~24 h incubation using the BD Phoenix™ M50 instrument and interpreted according to current Clinical and Laboratory Standards Institute. The results of ID and AST obtained using the conventional method were used as the standard for comparison.

### Rapid ID and AST using the short-term incubation method

After 5~7 h incubation, microbial ID was obtained directly from the growth of bacteria on BAP by Zybio EXS2000 MS. Colonies on the CAP were taken for AST detection when the rapid ID was shown as GN bacteria. However, for Gram positive (GP) bacteria, the colony on BAP was selected. The suspension obtained from bacterial growth for 5~7 h could be used for AST detection on BD Phoenix™ M50.

The comparisons between rapid and conventional methods were categorized as: category agreement (CA), essential agreement (EA), very major error (VME) (false susceptibility), major error (ME) (false resistance), or minor error (mE) (susceptible/resistant vs. intermediate).

In this study, 57 isolates were subjected to both rapid and conventional method for AST detection.

#### Evaluation of the turn-around time (TAT) in the conventional method

The turn-around time (TAT) consisted of two components: (i) time to positivity, corresponding to the time required for microorganism growth by Bactec incubation; and (ii) processing time, corresponding to the time required to generate the final report (including bacterial ID, AST, validation of results, and reporting to clinicians).

#### Quality control

Standard strains *Staphylococcus aureus* ATCC 29,213 and *Enterococcus faecalis* ATCC 29,212 were used for internal quality control of PMIC-92 plate. *Streptococcus pneumoniae* ATCC 49,619 was used as SMIC/ID-2 QC strains. In addition, standard *Escherichia coli* ATCC 25,922 and *Pseudomonas aeruginosa* ATCC 27,853 were used as internal quality control strains for NMIC-413.

## Results

#### Comparison of rapid and conventional ID methods

A total of 124 monomicrobial-positive BCs were collected, including 67 GN bacteria, 53 GP bacteria and 4 *Candida* spp. In addition, 5 samples that were polymicrobial after subculture were excluded. These monomicrobial-positive BCs were detected using the rapid ID method combined with MALDI-TOF MS and compared to conventional method. Among the 67 GN isolates, 61 (91.0%) showed a score higher than 2, while 5 (7.5%) scored between 1.7 and 2, 1 (1.5%) was unidentified

(Table 1). Among the 53 GP isolates, 38 (71.7%) demonstrated a score higher than 2, 10 (18.9%) demonstrated a score between 1.7 and 2, and 5 (9.4%) Coagulase-negative staphylococci (CoNS) were incorrectly identified as another CoNS (Table 2). Among 4 *Candida* spp, 3 (75%) demonstrated a score higher than 2, 1 (25%) demonstrated a score between 1.7 and 2 (Table 2).

Compared with the conventional method, the concordance of the rapid ID results of 67 GN bacteria was 98.5% (66/67), while the concordance of 53 GP bacteria was 90.6% (48/53) at the species level, but 98.1% (52/53) at the genus level. The *Candida* isolates exhibited perfect concordance rate of 100% (4/4).

#### Comparison of rapid and conventional AST results using BD M50 AST panels

Among the 124 isolates, 57 isolates were selected for both rapid and conventional AST using BD M50, including 37 *Enterobacteriaceae* (27 *E. coli*, 4 *K. pneumoniae*, 2 *C. freundii*, 2 *E. cloacae*, 1 *S. marcescens*, 1 *S. typhimurium*) and 5 non-fermenting gram-negative rods (1 *P. aeruginosa*, 1 *A. radioresistens*, 1 *A. pittii*, 1 *A. junii*, 1 *A. veronii*), 8 *staphylococcus*, 3 *Enterococcus*, 4 *Streptococcus*. 51.8% (14/27) of *E. coli* and 25.0% (1/4) of *K. pneumoniae* were resistant to Ceftriaxone, 2 isolates were Carbapenem-resistant *Enterobacteriaceae* (CRE), 3 isolates of *Staphylococci* were methicillin resistant staphylococcus (MRS). For all 42 GN isolates (except *A. junii*), a total of 946 bacterial-antimicrobial combinations were analyzed. Compared to the conventional AST, the EA, CA, mE, ME, VME of the BD M50 NMIC-413 panels by rapid AST were 98.5%, 97.1%, 2.3%, 0.4% and 0.1%, respectively

**Table 1** Gram-negative bacteria from monomicrobial blood cultures identified by rapid and conventional culture-dependent method ( $n=67$ )

Organisms	Conventional ID			Rapid ID				
	>2.0	1.7-2.0	<1.7	>2.0	1.7-2.0	<1.7	un-ID	mis-ID
<i>E. coli</i>	39	1		35	5			
<i>K. pneumoniae</i>	9			9				
<i>P. aeruginosa</i>	3			3				
<i>C. freundii</i>	2			2				
<i>S. marcescens</i>	2			2				
<i>E. cloacae</i>	2			2				
<i>A. junii</i>	1			1				
<i>B. cenocepacia</i>	1			1				
<i>A. pittii</i>	1			1				
<i>A. radioresistens</i>	1			1				
<i>P. mirabilis</i>	1			1				
<i>A. veronii</i>	1			1				
<i>A. xylosoxidans</i>	1			1				
<i>S. typhimurium</i>	1			1				
<i>Bru. spp</i>	1						1	
Total isolates	66	1		61	5		1	

un-ID: unidentified, mis-ID: misidentified

**Table 2** Gram-positive bacteria and fungi from monomicrobial blood cultures identified by rapid and conventional culture-dependent method ( $n=57$ )

Organisms	Conventional ID			Rapid ID			un-ID	mis-ID
	>2.0	1.7-2.0	<1.7	>2.0	1.7-2.0	<1.7		
<i>S. hominis</i>	17	1		11	5			2
<i>S. epidermidis</i>	7	1		5	2			1
<i>S. aureus</i>	6			6				
<i>S. capitis</i>	5			4	1			
<i>S. haemolyticus</i>	3			2				1
<i>S. warneri</i>	1				1			
<i>S. lugdunensis</i>	1	1			1		1	
<i>E. faecalis</i>	4			4				
<i>S. oralis</i>	1			1				
<i>S. sanguinis</i>	1			1				
<i>S. agalactiae</i>	1			1				
<i>S. dysgalactiae</i>	1			1				
<i>L. monocytogenes</i>	1			1				
<i>C. striatum</i>	1			1				
<i>C. glabrata</i>	2			1	1			
<i>C. tropicalis</i>	1			1				
<i>C. lusitanae</i>	1			1				
Total isolates	54	3		41	11		1	4

un-ID: unidentified, mis-ID: misidentified

**Table 3** AST Results obtained using the rapid method compared with those of the conventional method in Gram-negative bacteria (41) of NMIC-413

Antimicrobial agent	N	N (%) of				
		EA	CA	mE	ME	VME
Amikacin	41	41(100)	41(100)			
Gentamicin	41	41(100)	41(100)			
Tobramycin	41	39(95.1)	36(87.8)	4(9.8)	1(2.4)	
Ertapenem	37	37(100)	36(97.3)		1(2.7)	
Imipenem	41	39(95.1)	39(95.1)	2(4.9)		
Meropenem	41	39(95.1)	40(97.6)		1(2.4)	
Cefazolin	37	37(100)	36(97.30)	1(2.7)		
cefuroxime	37	37(100)	37(100)			
cefoxitin	37	37(100)	37(100)			
Ceftazidime	41	41(100)	41(100)			
Ceftriaxone	40	39(97.5)	40(100)			
Cefepime	41	40(97.6)	39(95.1)	2(4.9)		
Aztreonam	39	36(92.3)	36(92.3)	3(7.7)		
Amoxicillin/ clavulanic	37	36(97.3)	36(97.3)			1(2.7)
Ampicillin/sulbactam	39	39(100)	35(89.7)	4(10.3)		
Piperacillin/tazobactam	41	41(100)	39(95.1)	2(4.9)		
Sulfamethoxazole/trimethoprim	40	40(100)	40(100)			
chloramphenicol	40	39(97.5)	38(95)	1(2.5)	1(2.5)	
Ciprofloxacin	41	41(100)	41(100)			
Levofloxacin	41	40(97.6)	40(97.6)	1(2.4)		
Norfloxacin	37	37(100)	37(100)			
minocycline	39	39(100)	38(97.4)	1(2.6)		
tetracycline	40	40(100)	40(100)			
Tigecycline	37	37(100)	36(97.3)	1(2.7)		
Total	946	932(98.5)	919(97.1)	22(2.3)	4(0.4)	1(0.1)

**Table 4** AST Results obtained using the rapid method compared with those of the conventional method in *Staphylococcus* and *Enterococcus* of PMIC-92

Antimicrobial agent	N	N (%) of				
		EA	CA	mE	ME	VME
Gentamicin	8	8(100)	8(100)			
ceftaroline	4	4(100)	4(100)			
ampicillin	3	3(100)	3(100)			
penicillin	11	11(100)	11(100)			
Oxacillin	8	8(100)	8(100)			
daptomycin	11	11(100)	11(100)			
Sulfamethoxazole /trimethoprim	8	8(100)	8(100)			
teicoplanin	11	11(100)	11(100)			
Vancomycin	11	11(100)	11(100)			
Clindamycin	8	8(100)	8(100)			
Erythromycin	11	10(90.9)	10(90.9)	1(9.1)		
chloramphenicol	11	11(100)	10(90.9)	1(9.1)		
Linezolid	11	11(100)	11(100)			
Ciprofloxacin	11	11(100)	11(100)			
Levofloxacin	11	11(100)	11(100)			
Rifampin	8	8(100)	8(100)			
minocycline	11	11(100)	11(100)			
tetracycline	11	11(100)	11(100)			
Tigecycline	11	11(100)	11(100)			
Total	179	178(99.4)	177(98.9)	2(1.1)		

**Table 5** AST Results obtained using the rapid method compared with those of the conventional method in *Streptococci* of SMIC/ID-2

Antimicrobial agent	N	N (%) of				
		EA	CA	mE	ME	VME
Meropenem	4	4 (100)	4 (100)			
cefotaxime	4	4 (100)	4 (100)			
cefepime	4	4 (100)	4 (100)			
amoxicillin	4	4 (100)	4 (100)			
penicillin	4	4 (100)	4 (100)			
Vancomycin	4	4 (100)	4 (100)			
Clindamycin	4	4 (100)	4 (100)			
Erythromycin	4	4 (100)	4 (100)			
chloramphenicol	4	4 (100)	4 (100)			
Linezolid	4	4 (100)	4 (100)			
Levofloxacin	4	4 (100)	4 (100)			
tetracycline	4	4 (100)	4 (100)			
Total	48	48 (100)	48 (100)	0	0	0

(Table 3). For 8 *Staphylococcus* spp and 3 *Enterococcus* spp, a total of 179 bacterial-antimicrobial combinations were analyzed. The EA, CA, mE of PMIC-92 panels were 99.4%, 98.9% and 1.1%, respectively, while ME and VME were 0% (Table 4). For the SMIC/ID-2 panel, the EA, CA were both 100% of 4 *Streptococcus* spp (Table 5).

#### Evaluation of the TAT

During the study, we collected a total of 393 positive BC bottles, of which 264 (67.2%) were able to utilize the rapid ID and AST method. The mean TAT for detecting and reporting 57 isolates using conventional and rapid

method were  $73.0 \pm 10.7$  h and  $48.9 \pm 10.2$  h, respectively. The TAT of the rapid method was shortened by about 24 h (Table 6).

#### Discussion

Immediate administration of appropriate antibiotics is necessary for the effective treatment of bacteremia, as any delay is associated with increased morbidity and mortality [14]. MALDI-TOF MS has been proven to directly identify bacteria in positive BCs [2, 15, 16]. However, it has drawbacks such as time-consuming, complicated operation, and expensive [17, 18]. Molecular

**Table 6** Mean and standard deviation value of times in conventional method

Times (h)	Mean	GN	GP
Time to positivity	11.7 ± 4.8 (3.3 ~ 24.7)	10.7 ± 4.2 (3.3 ~ 18.7)	14.6 ± 5.2 (7.7 ~ 24.7)
<i>Conventional method</i>			
Processing time	61.3 ± 12.4 (27.4 ~ 87.6)	62.4 ± 11.8 (35.6 ~ 87.6)	57.9 ± 13.9 (27.4 ~ 79.5)
TAT	73.0 ± 10.7 (48.1 ~ 96.1)	73.1 ± 10.6 (48.3 ~ 96.1)	72.5 ± 11.5 (48.1 ~ 90.2)
<i>Rapid method</i>			
Processing time	37.2 ± 11.8 (14.1 ~ 65.2)	38.3 ± 11.7 (14.1 ~ 65.2)	34.2 ± 12.2 (16.0 ~ 57.1)
TAT	48.9 ± 10.2 (26.7 ~ 73.7)	49.0 ± 10.4 (26.7 ~ 73.7)	48.8 ± 9.8 (36.7 ~ 67.7)

Minimum and maximum values are reported in brackets

methods have shown to be efficient for the rapid ID of specific microorganisms but can identify a limited range of microorganisms and antimicrobial resistance genes. Nucleic acids-based technologies could be used for sensitive detection of bloodstream pathogens directly from a blood sample. Due to target limitations, most commercial platforms require a combination of standard blood cultures and adjunctive molecular detection. Additionally, all methods still require a culture step to obtain isolates for comprehensive AST [19]. Rapid AST methods have been developed in recent years, including EUCAST rapid AST and CLSI disk diffusion using positive blood culture medium [20–21]. However, the rapid AST methods require strain identification and cannot be performed in laboratories without mass spectrometry instruments. Moreover, their manual setup and the imperative requirement to read the inhibition zone diameters at strictly defined time points are tremendously labor-intensive [20]. Although the rapid AST has important value for some critically ill patients, its shortcomings hinder its large-scale use in clinical microbiology laboratories. In addition, the CLSI method detects a limited variety of bacterial strains, which cannot meet clinical needs.

The application of MALDI-TOF MS for identifying bacterial colonies from solid media has significantly improved and accelerated routine microbiological diagnostics [8, 22]. Some studies have described a new process based on the short-term incubation method for identifying bacterial pathogens, including ID and AST, from blood cultures [13, 23, 24]. One study showed that the optimal incubation time to ascertain GP bacterial ID was 4.5 h, but for GN bacterial ID was 3.5 h, the identification rates were 97.4% for GN bacteria and 100% for GP bacteria when compared to the conventional method [23]. Another study showed that the species-level ID concordance rate after 6 h of incubation was 90.9% (159/175), and 80.6% (141/175) after 4 h of incubation

[25]. To facilitate operation and standardization, we used the incubation time of 5~7 h for GP, GN bacteria and fungi. In this study, 9 bacterial strains were incubated for 5 h, 11 strains were incubated for 6 h, and 37 strains for 7 h, with an average incubation time of 6.5 h (median time of 7 h). There was no significant difference in the consistency of ID and AST between strains cultured for 5–7 h. Then, the rapid ID and AST of colonies were done by Zybion EXS2000 MS and BD M50, respectively. This means that processing positive BC samples, as well as ID and AST, can be completed on the same day.

With respect to rapid ID, our results demonstrate that the performance of the presented method is very high and satisfactory for both GP and GN isolates. In our study, only 1 GN bacteria was unidentified, it turned out to be *Brucella* spp with a time to positive of nearly 4 days. The cause of this discrepancy might be the very slow growth of *Brucella*, which could result in a failure to produce sufficient proteins for accurate MALDI-TOF MS analysis. Besides, 1 strain of *Enterobacter cloacae* was misidentified as *Enterobacter asburiae*. Since both belong to the *Enterobacter cloacae* complex, we consider that the ID results are consistent. Low identification scores (<2.0) were mainly associated with CoNS such as *Staphylococcus epidermidis*, *Staphylococcus hominis*, and *Staphylococcus capitis*. A total of 5 GP isolates were incorrectly identified, including 1 strain of *Staphylococcus lugdunensis* was unidentified, and 4 CoNS strains were misidentified (including 2 *S. hominis*, 1 *S. haemolyticus* and 1 *S. epidermidis*), but they were consistent at genus level. No CoNS were mis-identified as *S. lugdunensis*, which is virulent and the symptoms are similar to that of *S. aureus*. However, CoNS were usually recognized as contaminant species in blood samples, our method showed a greater accuracy when excluding these contaminant strains. The concordance rate of fungi ID was 100%. Moreover, previous studies have shown that the ID rate of GN bacteria is higher than that of GP bacteria [26, 27], and our study has reached similar conclusions.

Regarding the applications of this rapid procedure for AST, positive BCs for GN bacteria were assessed with 24 antimicrobial agents using NMIC-413. The rates of EA, CA, mE, ME, VME were 98.5%, 97.1%, 2.3%, 0.4% and 0.1%, respectively. Meanwhile, 11 isolates of GP were evaluated for 19 antimicrobial agents using PMIC-92, the rates of EA, CA, mE were 99.4%, 98.9% and 1.1%. Similar findings have been reported in previous studies [13, 24]. It is noteworthy that 4 isolates of *Streptococci* were assessed for 12 antimicrobial agents using SMIC/ID-2, both EA and CA exhibited flawless outcomes with a rate reaching 100%. The rapid AST showed excellent results with low error percentage (mE, ME, VME) meeting the performance standards of the U.S. Food and Drug Administration (FDA) (mE<3%, VME<1.5%) [28]. It

should be noted that colonies on the CAP were used for rapid AST in GN bacteria to avoid contamination from other GP bacteria. One *Acinetobacter junii* was excluded due to the AST without acquired MIC in both rapid and conventional methods, and finally results of AST were obtained by KB method. ME and VME were found only in 3 isolates (1 *Escherichia coli*, 1 *Enterobacter cloacae*, 1 *Salmonella typhimurium*). On the other hand, mE was mostly observed with Ampicillin/sulbactam, Tobramycin, Aztreonam. In the isolated *Staphylococci*, 100% EA and CA were detected for oxacillin, indicating that our method is suitable for detecting methicillin-resistant *Staphylococcus* (MRS). Finally, in addition to high performance rates for ID and AST, the rapid method greatly shortens the processing time, resulting in the reporting time being about 24 h earlier than traditional methods (shortening TAT by nearly one-third). The results of rapid method have been recognized by many clinicians in our hospital.

The rapid method used in this study also has some limitations. The method cannot be used for polymicrobial-positive BCs. Approximately 3.9% (5/129) of the positive BCs collected in this study were detected as polymicrobial, which is slightly lower than other studies [24, 25]. Besides, the commercialized antimicrobial susceptibility testing panel used in our laboratory do not contain novel  $\beta$ -lactam combination agents, such as imipenem-relebactam, meropenem-vaborbactam, Ceftazidime-avibactam, and no AST was performed on these new drugs in this study. However, there are several highlights in our article. Firstly, the short-term culture 5~7 h on solid medium method for positive BCs showed a high concordance (>90%). Secondly, the rapid method exhibited excellent performance in AST using different AST panels of BD-M50 system. Thirdly, this new procedure allows for a reduction in TAT by nearly 24 h. In summary, the rapid method can be applied to most positive BC samples and provides a workflow friendly approach that can obtain results faster and more reliably, which is very beneficial for the treatment of sepsis patients.

#### Acknowledgements

Not applicable.

#### Author contributions

PPT, WF and HWY designed and wrote the manuscript. SSS, LSZ and TW was responsible for obtaining data in the hospital's HIS and LIS systems. PPT, HY, MYD, CZD, LW and HWY collected the data and completed the analysis. PPT, WF and HWY participated in the manuscript preparation. PPT and HWY participated in the study design and data analysis. All authors reviewed the results and data analysis and contributed comments. All authors agreed on the final version for submission to the journal.

#### Funding

This work was supported by the Natural Science Foundation of Hubei Province, China [Grant Number 2021CFB261 (to H-WY)], Jingzhou science and technology plan projects, Hubei Province, China [grant number 2023HC82 (to

H-WY)] and the Doctoral Research Start-up Fund Project of the First People's Hospital of Jingzhou, China [Grant Number 2023DIF03 (to H-WY)].

#### Data availability

The datasets used and analyzed during the current study are available from the corresponding author on reasonable request.

#### Declarations

##### Ethics approval and consent to participate

This study involving humans was approved by the Ethics Committee of the First Affiliated Hospital of Yangtze University, Jing Zhou, China (no. LL2023124). Written informed consent was obtained from individual or guardian participants.

##### Consent for publication

Not applicable.

##### Competing interests

The authors declare no competing interests.

Received: 6 December 2023 / Accepted: 4 June 2024

Published online: 06 June 2024

#### References

- Adhikari NK, Fowler RA, Bhagwanjee S, Rubenfeld GD. Critical care and the global burden of critical illness in adults. *Lancet*. 2010;376(9749):1339–46.
- Morgenthaler NG, Kostrzewa M. Rapid identification of pathogens in positive blood culture of patients with sepsis: review and meta-analysis of the performance of the sepsityper kit. *Int J Microbiol*. 2015;2015:827416.
- Goto M, Al-Hasan MN. Overall burden of bloodstream infection and nosocomial bloodstream infection in North America and Europe. *Clin Microbiol Infect*. 2013;19(6):501–9.
- Seifert H. The clinical importance of microbiological findings in the diagnosis and management of bloodstream infections. *Clin Infect Dis*. 2009;48:S238–45.
- Koupetori M, Retsas T, Antonakos N, Vlachogiannis G, Perdios I, Nathanail C, et al. Bloodstream infections and sepsis in Greece: over-time change of epidemiology and impact of de-escalation on final outcome. *BMC Infect Dis*. 2014;14:272–81.
- Kumar A, Roberts D, Wood KE, Light B, Parrillo JE, Sharma S, et al. Duration of hypotension before initiation of effective antimicrobial therapy is the critical determinant of survival in human septic shock. *Crit Care Med*. 2006;34:1589–96.
- Kerremans JJ, Verboom P, Stijnen T, Hakkaart-van Roijen L, Goessens W, Verbrugh HA, et al. Rapid identification and antimicrobial susceptibility testing reduce antibiotic use and accelerate pathogen-directed antibiotic use. *J Antimicrob Chemother*. 2008;61(2):428–35.
- Vlek AL, Bonten MJ, Boel CH. Direct matrix-assisted laser desorption ionization time-of-flight mass spectrometry improves appropriateness of antibiotic treatment of bacteremia. *PLoS ONE*. 2012;7:e32589.
- van Belkum A, Chatellier S, Girard V, Pincus D, Deol P, Dunne WM. Progress in proteomics for clinical microbiology: MALDI-TOF MS for microbial species identification and more. *Expert Rev Proteom*. 2015;12(6):595–605.
- de Almeida JN, Sztajn bok J, da Silva AR, Vieira VA, Galastrì AL, Bissoli L, et al. Rapid identification of moulds and arthroconidial yeasts from positive blood cultures by MALDI-TOF mass spectrometry. *Med Mycol*. 2016;54:885–9.
- Faron ML, Buchan BW, Ledebøer NA. Matrix-assisted laser desorption ionization-time of Flight Mass Spectrometry for Use with positive blood cultures: methodology, performance, and optimization. *J Clin Microbiol*. 2017;55(12):3328–38.
- Costa A, Figueroa-Espinosa R, Martínez JA, Fernández-Canigia L, Maldonado MI, Bergese SA, et al. MALDI-TOF MS-Based KPC Direct Detection from patients' positive blood culture bottles, short-term cultures, and colonies at the hospital. *Pathogens*. 2023;12(7):865–83.
- Ha J, Hong SK, Han GH, Kim M, Yong D, Lee K. Same-day identification and Antimicrobial susceptibility testing of Bacteria in positive blood culture broths using short-term incubation on solid medium with the MicroFlex LT, Vitek-MS, and Vitek2 systems. *Ann Lab Med*. 2018;38:235–41.

14. Ascione T, Pagliano P, Mariconda M, Rotondo R, Balato G, Toro A, et al. Factors related to outcome of early and delayed prosthetic joint infections. *J Infect*. 2015;70:30–6.
15. Schneiderhan W, Grundt A, Wörner S, Findeisen P, Neumaier M. Work flow analysis of around-the-clock processing of blood culture samples and integrated MALDI-TOF mass spectrometry analysis for the diagnosis of blood-stream infections. *Clin Chem*. 2013;59:1649–56.
16. Pan HW, Li W, Li RG, Li Y, Zhang Y, Sun EH. Simple sample preparation method for direct microbial identification and susceptibility testing from positive blood cultures. *Front Microbiol*. 2018;9:481.
17. Monteiro J, Inoue FM, Lobo AP, Sugawara EK, Boaretti FM, Tufik S. Fast and reliable bacterial identification direct from positive blood culture using a new TFA sample preparation protocol and the Vitek® MS system. *J Microbiol Methods*. 2015;109:157–9.
18. Jo SJ, Park KG, Han K, Park DJ, Park YJ. Direct identification and antimicrobial susceptibility testing of bacteria from positive blood culture bottles by matrix-assisted laser desorption/ionization time-of-flight mass spectrometry and the Vitek2 system. *Ann Lab Med*. 2016;36:117–23.
19. Sinha M, Jupe J, Mack H, Coleman TP, Lawrence SM, Fraley SI. Emerging technologies for molecular diagnosis of sepsis. *Clin Microbiol Rev*. 2018;31(2):e00089–17.
20. The European Committee on Antimicrobial Susceptibility Testing. Zone diameter breakpoint tables for rapid antimicrobial susceptibility testing (RAST) directly from blood culture bottles. Version 6.1, 2023. <http://www.eucast.org>.
21. Clinical and Laboratory Standards Institute. CLSI M02. Performance Standards for Antimicrobial Disk Susceptibility Tests. 14th Edition. Wayne, PA: CLSI, 2024.
22. Huang AM, Newton D, Kunapuli A, Gandhi TN, Washer LL, Isip J, et al. Impact of rapid organism identification via matrix-assisted laser desorption/ionization time-of-flight combined with antimicrobial stewardship team intervention in adult patients with bacteremia and candidemia. *Clin Infect Dis*. 2013;57:1237–45.
23. Tsai YW, Lin TC, Chou HY, Hung HY, Tan CK, Wu LC, et al. Shortening the time of the identification and Antimicrobial susceptibility testing on positive blood cultures with MALDI-TOF MS. *Diagnostics (Basel)*. 2021;11:undefined.
24. Mauri C, Principe L, Bracco S, Meroni E, Corbo N, Pini B, et al. Identification by mass spectrometry and automated susceptibility testing from positive bottles: a simple, rapid, and standardized approach to reduce the turnaround time in the management of blood cultures. *BMC Infect Dis*. 2017;17:749.
25. Hong SK, Chang BK, Song SH, Kim EC. Use of MALDI-TOF MS technique for rapid identification of bacteria from positive blood cultures. *Indian J Med Microbiol*. 2014;32:419–22.
26. Idelevich EA, Schüle I, Grünastel B, Wüllenweber J, Peters G, Becker K. Rapid identification of microorganisms from positive blood cultures by MALDI-TOF mass spectrometry subsequent to very short-term incubation on solid medium. *Clin Microbiol Infect*. 2014;20:1001–6.
27. Altun O, Botero-Kleiven S, Carlsson S, Ullberg M, Özenci V. Rapid identification of bacteria from positive blood culture bottles by MALDI-TOF MS following short-term incubation on solid media. *J Med Microbiol*. 2015;64:1346–52.
28. Center for Devices and Radiological Health, Food and Drug Administration. Class II Special Controls Guidance Document: Antimicrobial Susceptibility Test (AST) Systems, <https://www.fda.gov/downloads/MedicalDevices/Device-RegulationandGuidance/GuidanceDocuments/ucm388961.pdf> (Updated on Aug 2009).

### Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.